




GENICON[®]

Suction Irrigation




Sturdy probe for organ manipulation quickly connects with handpiece. Handpiece also fits many competitive probes.

Grip feature built into handpiece

Rear port for easy access inside chamber

Product and tubing are DEHP free



Gravity Fed Disposable Suction / Irrigation Tube Sets

| Description | Size | Qty/Box | Part Number |
|--|------------|---------|-------------|
| 8-Vent Disposable Probe with Gravity Fed Tube Set (Single pouch, bag spike) | 5mm x 33cm | 10 | 700-000-010 |

Powered Disposable Suction / Irrigation Tube Sets

| Description | Size | Qty/Box | Part Number |
|--|------------|---------|--------------|
| 8-Vent Disposable Probe and Powered Tube Set (Single pouch, bag spike, and clamp) | 5mm x 33cm | 10 | 710-300-000* |

*Only for Customers with Discontinued Pump 717-000-000

No External Pump or Battery


Unique one-piece design separates sealed battery from saline bag to prevent electrical shorting

Color Coded Buttons

Ergonomic design for ambidextrous toggling of suction and irrigation functions

Hand-piece Compatibility

Standard thread pattern works with a variety of disposable and reusable probes*



Powerful Flow Rate

Variable speed dynamic flow provides up to 4L/min of focused lavage to clear the operative site and surrounding tissue structures

Multi-Grip Handle

User-friendly design provides comfort and control at many grip positions

Product and tubing are DEHP free



Side Vents

Side vents on standard probes allow for strong, controlled lateral spray

Battery Powered Disposable Suction / Irrigation Tube Sets

| Description | Size | Qty/Box | Part Number |
|--|---------------------------|---------|-------------|
| Aquas PowerFlow Tube Set (No Probe) | --- | 8 | 720-005-000 |
| Aquas PowerFlow Tube Set (One Probe) | 5mm x 33cm | 8 | 720-005-001 |
| Aquas PowerFlow Tube Set (Two Probes) | 5mm x 33cm 10mm x 33cm | 8 | 720-005-002 |

*It is the responsibility of the end user to verify all instrument compatibility from other manufacturers prior to performing any procedure. Follow all instructions for use.

Reusable Vent Probes

| Description | 5mm x 34cm | 5mm x 45cm | 10mm x 34cm |
|---------------|-------------|-------------|-------------|
| 16-Vent Probe | 707-005-033 | 707-005-045 | 707-010-033 |

Reusable Electrocautery Probes

| Description | Spatula Probe | J-Hook Probe | L-Hook Probe |
|-----------------------------------|---------------|--------------|--------------|
| Electrocautery Probe (5mm x 33cm) | 707-005-133 | 707-005-233 | 707-005-333 |



Olga Haberland

International Business Development / Marketing
 Phone: +1 407.657.4851 x304
 Email: olh@geniconendo.com
 Skype: ohaberland
 Languages: English, Spanish

Refaat Hafez, MD

Middle East and Africa Market Manager
 Phone: +1 407.657.4851 x303
 Fax: +1 407.679.0880
 Email: rmonem@geniconendo.com
 Skype: refaatmoniem
 Languages: English, Arabic

Justin Colon

Customer Service
 Phone: +1.407.657.4851 x800
 E-mail: customerservice@geniconendo.com
 Skype: jcc_450

Laura Strelitzki

Latin America and Asia Business Development
 Phone: +1.407.657.4851 x308
 E-mail: lal@geniconendo.com
 Skype: laura.acosta.lopez1
 Languages: English, Spanish

Doug Dotson

Clinical Product Specialist
 Phone: +1 407.848.9728
 Secondary Phone: +1.407.657.4851 x800
 Email: dwd@geniconendo.com
 Skype: dwd_102

Jim Kirchberg

Vice President of Sales - USA
 Phone: +1 815.355.5603
 Secondary Phone: +1.407.657.4851 x302
 Email: JimK@geniconendo.com
 Skype: jim.kirchberg.genicon

Visit www.geniconendo.com for a complete list of trademarks and patents