



GENICON[®]


Specimen Retrieval Systems




All-in-One Introducer Style Specimen Bags

Laparoscopic Solutions: Specimen Retrieval Bag

The GeniStrong™ Bag is an easy-to-use retrieval system designed to temporarily contain specimen and facilitate its removal from the patient during laparoscopic surgery. Designed to minimize contamination of the abdominal cavity, the Specimen Retrieval Bag includes the following distinctive features:



Burst Strength Comparison



*Test Report Data On File

Actual Size

550-000-008

2000mL
14cm x 27cm
15mm Introducer

550-000-005

750mL
11cm x 22cm
12mm Introducer

550-000-003

240mL
7.3cm x 17cm
10mm Introducer

GENICON

GENICON

GENICON

**Play to Your Strengths
Ask for the Strong Bag**



Bag-Only Specimen Bags



Length	Mouth Opening	Qty/Box	Volume	Part Number
3 in. / 7 cm	3.5 in. / 9 cm	20	50 mL	550-000-017
5 in. / 13 cm	3.5 in. / 9 cm	20	150 mL	550-000-009
6.5 in. / 17 cm	4 in. / 11 cm	20	240 mL	550-000-010
8 in. / 20 cm	8 in. / 20 cm	10	750 mL	550-000-011
10 in. / 23 cm	9 in. / 25 cm	10	1600 mL	550-000-012
13 in. / 34 cm	10 in. / 25 cm	10	3000 mL	550-000-013

Features & Benefits

Feature	GENICON	Anchor	Applied	Medtronic	Ethicon
Rip-Stop Nylon Bag Material	✓	✓		✓	
Superior Strength	✓	✓			
Will Not Burst / Propagate Tears	✓	✓		✓	
Multi-Use (Single Patient)	✓	✓			
Available on U.S. GPO Contracts	✓		✓	✓	✓
Cost Efficient	✓		✓		
Multiple Sizes Available	✓	✓	✓	✓	
Gas Seal in Introducer	✓		✓		
Buttonhole Quick Bag Release Feature	✓				
Metal Pull/Push Rod (Sturdy Placement)	✓	✓	✓		✓
Radiopaque Feature Visible on X-Rays	✓				
Tapered Bag Shape for Easy Removal	✓	✓	✓	✓	✓
Translucent	✓	✓	✓	✓	✓

Superior Strength

Superior Quality Provides Superior Performance
Does Not Burst Or Propagate Tears Like Poly (Plastic) Bags
Ideal For MIS Procedures
(Removing Large Specimen From Small Incisions)

Clinically Preferred

Improved Patient Outcomes
Improved Patient Safety
Industry 1st Radiopaque Feature Allows Detection On X-rays
Improved Patient And Surgeon Satisfaction

Rip-stop Nylon

Rip-stop Nylon With Polymer Coating For Leak-proof Performance
Improved Quality Of Care
Superior Strength Vs. Poly (Plastic) Bags

Cost Efficient


Available On Multiple GPO Contracts
Competitively Priced For The Health Care Provider
Reduced Procedural Costs
Can Be Used Multiple Times On A Single Patient

Recommended Anatomy Removal

Appendix, Ectopic Pregnancy, Gallbladder (Standard and Acute), Large Fibroid, Node Biopsy, Ovary, Prostate, Small Bowel Resection, Uterine Fibroids, VATS	550-000-003	550-000-010	550-000-000
Bowel Resection, Gastric Sleeve, Kidney, Lobe Resection (Liver/Lung), Spleen, Uterus, VATS	550-000-005	550-000-011	
Enlarged Spleen, Large Lobe Resection, Large Uterus, Malignant Kidney, Hemi-Colectomy, VATS	550-000-008	550-000-012 550-000-013	
Nanosurgery, Lymph Node Dissection, Oophorectomy, Pediatric Procedures	N/A	550-000-017 550-000-009	



Poly Introducer Style Specimen Bags




550-000-000

230mL

5.1cm x 12cm

10mm Introducer

Burst Strength Comparison



*Test Report Data On File