



HYPERION[®]

Large UVC Chamber



Rapid UVC disinfection at the touch of a button

Every Surface Exposed

Six plane exposure of high intensity UVC light

Easy to Use

No heat, no chemicals, no consumables

Rapid Disinfection

30 Second disinfection cycle



MOBILE
CAPABILITIES

GREATER
CAPACITY

LOWER COSTS

The leading design for today's demanding
total room disinfection.

www.finsentech.com

 **Finsentech**

UNIQUE DESIGN

The HYPERION® chamber is a unique surface decontaminator, effective against all known problematic pathogens. Its unique, patent-protected six plane ultraviolet (UVC) decontamination process penetrates hard-to-reach areas of products that are often missed when disinfecting by hand.

GUARANTEED SAFETY

The HYPERION® chamber is a full enclosure device, so there is no risk of operator exposure to UVC. As an added safety feature, a digitally controlled TITAN interlock is fully integrated into the door mechanism.

SIMPLE TO USE

Simply open the door, place your products inside, close the door, and press 'GO'. The system will then perform a 30-second disinfection cycle.

EASY TO DEPLOY

Due to the HYPERION® chamber using only UVC light to disinfect, there is no requirement for specialist ventilation or extraction.

PLUG AND PLAY

HYPERION® only requires a standard single-phase mains power connection, so can be used in any location.

BUILT TO LAST

Constructed of aircraft-grade stainless steel, and with a modular design to allow easy access for routine maintenance, HYPERION® is designed to last indefinitely and does not have a fixed life span.

PROVEN TECHNOLOGY

Using UVC light at a wavelength of 253.7nm, across all six planes and at a known level of intensity/distance, ensures that any product placed into the HYPERION® chamber will be terminally disinfected.

MORE VALUE

HYPERION® does not require any chemicals of any kind, and there are no ongoing consumable costs.

HIGH THROUGHPUT

As the HYPERION® chamber only requires a 30-second cycle to terminally disinfect a product, it is capable of supporting a high capacity of disinfections within one single working day.

ADAPTABLE

HYPERION® has been engineered to be as adaptable as possible, and the chamber can take either multiple products or extremely large products in one single cycle due to the internal dimensions of 975mm x 650mm x 600mm.

If you continue to be concerned about surface and airborne pathogens, have been waiting for the next generation of UVC disinfection technology or if you have already acquired a UV system, HYPERION® can introduce you to entirely new levels of disinfection. Efficacy, Speed, Confidence, and Value. To learn more about how HYPERION® can help you make a better environment for your patients and professionals, please contact us at:

SALES@FinsenTech.com or 0207 193 5641



Fast

60-second cycle creates major time savings and facilitates a rapid turnaround of equipment.



Proven

Proven to kill pathogens.



Safe

No consumables needed and no harmful chemicals or residues.



Green

Works on sensitive equipment at room temperature and pressure with no risk of moisture ingress.



Validated

Provides a standardised, high-level decontamination solution.



Cost Effective

No need for manual cleaning.