

KEBOMED



Enteral Division

AMT Enteral Devices,
Bridles Accessories,
Surgical Devices and
Bladder & Bowel
Management
Devices

TABLE OF CONTENTS

AMT Bridle Nasal Tube Retaining System Family	3
Bridle Pro® AMT Bridle	
Mini Classic Gastrostomy Tube	4
MiniONE® Gastrostomy Tube Family	5-6
MiniONE® Balloon Button	
MiniONE® Capsule Non-Balloon Button	
Traditional Length Gastrostomy Tube Family	7-8
Balloon G-Tube	
Capsule Monarch® G-Tube (non-balloon)	
G-JET® Gastric-Jejunal Feeding Tube Family	9-13
micro G-JET® Button	
G-JET® Button	
Traditional G-JET®	
Accessories	14-16
Extension Sets	
CINCH® Tube Securement Device	
AMT Clamp Tube Securement Device Guidewires	
Transition Adapters: ENFit®	
Replacement Y-Port Adapters	
Surgical Devices	17
Initial Placement Gastrostomy Kits	
Bladder & Bowel Management	18-19
Devices MiniSPC and Drainage	
sets MiniACE® Enema Button	
Contact Us and References	20

MEDICAL DEVICE SOLUTIONS:

Applied Medical Technology, Inc. (AMT) is a global leader in the development, design, and manufacture of enteral feeding devices, and accessories.

We continue to introduce innovative products and materials, designed to increase patient comfort while improving patient outcomes. In our state-of-the-art facility near Cleveland, Ohio, we believe in manufacturing the highest quality medical devices to meet the needs of our customers and patients.

AMT has been a pioneering force within the medical device market, revolutionizing the standard of care for both pediatric and adult patients. Our Founder/CEO, Dr. Picha, co-invented and manufactured the first low profile gastrostomy button, changing the enteral nutrition community forever.

Over time, AMT's product portfolio has evolved to include nasal tube retention systems, innovative solutions for paediatric jejunal feeding, safe tube securement options, bowel management, and a complete line of ENFit® products.



Bridle Pro®

AMT Bridle Nasal Tube Retaining System Family



An AMT Exclusive. The New Standard of Care.

Flexible Tubing

Easy to clean and comfortable for patients.

Alignment Markings

Assures device alignment for proper placement.

Pre-Attached Clip

Improved patient safety and improved handling.

Ease of Placement

Magnets at the distal end help facilitate a quick and easy connection.

Bridle Family Clips

Clips accommodate several nasal tube sizes.



AMT BRIDLE PRO® PAEDIATRIC SIZE NASAL TUBE RETAINING SYSTEM

DESCRIPTION	ORDER CODE	BOX/QTY
5F-6F Paediatric Bridle Pro® Range (Purple)	4-420506P	5
8F-10F Paediatric Bridle Pro® Range (Teal)	4-420810P	5
8F Paediatric Bridle Pro® (White)	4-4208P	5

Kit contents include: flexible bridle catheter, stylet guide, retrieval probe, pre-attached clip, clip opening tool and water soluble lubricant.

AMT BRIDLE PRO® NASAL TUBE RETAINING SYSTEM

8F-10F Bridle Pro® Range (Teal)	4-420810	5
12F-14F Bridle Pro® Range (Dark Blue)	4-421214	5
16F-18F Bridle Pro® Range (Yellow)	4-421618	5
8F Bridle Pro® (White)	4-4208	5
10F Bridle Pro® (Teal)	4-4210	5
12F Bridle Pro® (Dark Blue)	4-4212	5


Kit contents include: flexible bridle catheter, stylet guide, retrieval probe, pre-attached clip, clip opening tool and water soluble lubricant.

AMT BRIDLE PRO® SPARE CLIPS

5F-6F Bridle Pro® Range (Purple)	4-460506	5
8F-10F Bridle Pro® Range (Teal)	4-460810	5
12F-14F Bridle Pro® Range (Dark Blue)	4-461214	5
16F-18F Bridle Pro® Range (Yellow)	4-461618	5
8F Bridle Pro® (White)	4-4608	5
10F Bridle Pro® (Teal)	4-4610	5
12F Bridle Pro® (Dark Blue)	4-4612	5

MiniClassic® Balloon Button


Low Profile Balloon Button Gastrostomy Tube



Reduce Leaks
Improve Comfort
Balloon Durability

ENFIT® ORDER #						
LENGTH (cm)	12F	14F	16F	18F	20F	24F
0.8	5-1208-I	5-1408-I				
1	5-1210-I	5-1410-I	5-1610-I	5-1810-I	5-2010-I	
1.2	5-1212-I	5-1412-I	5-1612-I	5-1812-I	5-2012-I	
1.5	5-1215-I	5-1415-I	5-1615-I	5-1815-I	5-2015-I	5-2415-I
1.7	5-1217-I	5-1417-I	5-1617-I	5-1817-I	5-2017-I	5-2417-I
2	5-1220-I	5-1420-I	5-1620-I	5-1820-I	5-2020-I	5-2420-I
2.3	5-1223-I	5-1423-I	5-1623-I	5-1823-I	5-2023-I	
2.5	5-1225-I	5-1425-I	5-1625-I	5-1825-I	5-2025-I	5-2425-I
2.7	5-1227-I	5-1427-I	5-1627-I	5-1827-I	5-2027-I	5-2427-I
3	5-1230-I	5-1430-I	5-1630-I	5-1830-I	5-2030-I	5-2430-I
3.5	5-1235-I	5-1435-I	5-1635-I	5-1835-I	5-2035-I	5-2435-I
4	5-1240-I	5-1440-I	5-1640-I	5-1840-I	5-2040-I	5-2440-I
4.4		5-1444-I	5-1644-I	5-1844-I	5-2044-I	5-2444-I
5		5-1450-I				
5.5		5-1455-I				
6		5-1460-I		5-1860-I		
6.5		5-1465-I				

MiniClassic® Button comes packaged as 1/box (sterile). Kit contents include: device (1), bolus extension set (1), y-port extension set (1), 5cc/10cc balloon syringe (1), 35cc bolus syringe (1)




MiniONE® Balloon Button

Low Profile Balloon Button Gastrostomy Tube

Study Confirms: Over 70% of families prefer the AMT MiniONE® over the Avanos Mic-Key®.⁴

Balloon Vol. Clearly Labeled

The Balloon fill volume printed directly on balloon fill port.




Feed Sets Included in Kit

MiniONE® Button comes packaged as 1/box (sterile)
Kit contents includes: Device, bolus extension set,
Y-port extension set, 5ml, 10ml syringe, introducer
(12-18F) gauze, water soluble lubricant.

MiniONE® Balloon Button | Order Information

ENFit® ORDER #						
LENGTH (cm)	12F	14F	16F	18F	20F	24F
0.8	M1-5-1208-I	M1-5-1408-I				
1	M1-5-1210-I	M1-5-1410-I	M1-5-1610-I	M1-5-1810-I	M1-5-2010-I	
1.2	M1-5-1212-I	M1-5-1412-I	M1-5-1612-I	M1-5-1812-I	M1-5-2012-I	
1.5	M1-5-1215-I	M1-5-1415-I	M1-5-1615-I	M1-5-1815-I	M1-5-2015-I	M1-5-2415-I
1.7	M1-5-1217-I	M1-5-1417-I	M1-5-1617-I	M1-5-1817-I	M1-5-2017-I	M1-5-2417-I
2	M1-5-1220-I	M1-5-1420-I	M1-5-1620-I	M1-5-1820-I	M1-5-2020-I	M1-5-2420-I
2.3	M1-5-1223-I	M1-5-1423-I	M1-5-1623-I	M1-5-1823-I	M1-5-2023-I	M1-5-2423-I
2.5	M1-5-1225-I	M1-5-1425-I	M1-5-1625-I	M1-5-1825-I	M1-5-2025-I	M1-5-2425-I
2.7	M1-5-1227-I	M1-5-1427-I	M1-5-1627-I	M1-5-1827-I	M1-5-2027-I	M1-5-2427-I
3	M1-5-1230-I	M1-5-1430-I	M1-5-1630-I	M1-5-1830-I	M1-5-2030-I	M1-5-2430-I
3.5	M1-5-1235-I	M1-5-1435-I	M1-5-1635-I	M1-5-1835-I	M1-5-2035-I	M1-5-2435-I
4	M1-5-1240-I	M1-5-1440-I	M1-5-1640-I	M1-5-1840-I	M1-5-2040-I	M1-5-2440-I
4.4		M1-5-1444-I	M1-5-1644-I	M1-5-1844-I	M1-5-2044-I	M1-5-2444-I
5		M1-5-1450-I	M1-5-1650-I	M1-5-1850-I	M1-5-2050-I	M1-5-2450-I
5.5		M1-5-1455-I	M1-5-1655-I	M1-5-1855-I	M1-5-2055-I	M1-5-2455-I
6		M1-5-1460-I	M1-5-1660-I	M1-5-1860-I	M1-5-2060-I	M1-5-2460-I
6.5		M1-5-1465-I				
7		M1-5-1470-I				
8		M1-5-1480-I				
9		M1-5-1490-I				
10		M1-5-14100-I				



MiniONE® Capsule Non-Balloon Button

Low Profile Non-Balloon Gastrostomy Tube



Exclusive Capsule Technology

Material facilitates easy and comfortable insertion.

Glow Green Technology

Easier nighttime feeds.

Improve Comfort

Combines the advantage of a low profile device with the longer life of a non-balloon device.

Exclusive AMT Capsule Technology




More Durable

Internal bolster is smaller than most devices, durable and less likely to block the pylorus.¹²

Feed Sets Included in Kit

MiniONE® Button comes packaged as 1/box (sterile). Kit contents include: device (1), bolus extension set (1), y-port extension set (1), 10cc syringe (1), 35cc bolus syringe (1), Placement T-Handle/ Removal T-Handle‡, snap arm assembly, and gauze.



MiniONE® Balloon Button | Order Information

	14F	18F	20F	24F
LENGTH (cm)	ENFIT® ORDER #	ENFIT® ORDER #	ENFIT® ORDER #	ENFIT® ORDER #
1	M1-2-1410-I			
1.2	M1-2-1412-I	M1-2-1812-I		
1.5	M1-2-1415-I	M1-2-1815-I		
1.7	M1-2-1417-I	M1-2-1817-I	M1-2-2017-I	M1-2-2417-I
2	M1-2-1420-I	M1-2-1820-I		
2.4		M1-2-1824-I	M1-2-2024-I	M1-2-2424-I
2.5	M1-2-1425-I			
3	M1-2-1430-I	M1-2-1830-I		
3.4		M1-2-1834-I	M1-2-2034-I	M1-2-2434-I
4.4		M1-2-1844-I	M1-2-2044-I	M1-2-2444-I

Spare Removal Tools | For Obturated Non-Balloon and Capsule Non-Balloon Buttons

AMT Obturator for Distending button during insertion and removal.




‡ Non-Balloon Kit: Placement and Removal T-Handle

Capsule Non-Balloon Kit: Removal T – 6 Handle

FRENCH SIZE	DESCRIPTION	ORDER #	BOX/ QTY
14F	1.0, 1.2 cm	R-S1012	1
14F	1.5, 1.7 cm	R-S1517	1
14F	2.0 cm	R-S2023	1
14F	2.5 cm	R-S2427	1
14F	3.0 cm	R-S3033	1
18F	1.2 cm	R-L1012	1
18F, 20F, 24F	1.5, 1.7 cm	R-L1517	1
18F	2.0 cm	R-L2023	1
18F	2.4 cm	R-L2427	1
18F	3.0 cm	R-L3033	1
18F, 20F, 24F	3.4 cm	R-L3437	1
18F, 20F, 24F	4.4 cm	R-L4044	1

Traditional Balloon G-Tube

Adjustable Length Gastrostomy Tube



Adjustable External Bolster

Can accommodate various stoma lengths and has stabilizer feet that allow for air flow.

DESCRIPTION	ENFit® ORDER #
12F Traditional Balloon Gastrostomy Tube	7-1210-I
14F Traditional Balloon Gastrostomy Tube	7-1410-I
16F Traditional Balloon Gastrostomy Tube	7-1610-I
18F Traditional Balloon Gastrostomy Tube	7-1810-I
20F Traditional Balloon Gastrostomy Tube	7-2010-I
24F Traditional Balloon Gastrostomy Tube	7-2410-I

Reduce Leaks

- Exclusive AMT Balloon expands horizontally to provide superior stoma coverage.⁷

Improve Comfort + Safety

- Balloon protects the distal tip of the tube, lowering the risk for gastric ulceration.
- Enhanced External Retention Ring minimizes contact and allows skin to breathe.

Ease of Placement

- Often placed in the office or on-site.
- Tapered tip design for easier placement.
- Guidewire compatible.
- Radiopaque for x-ray confirmation.



AMT Balloon G-Tube is packaged 1/Box. All ENFit® Balloon G-Tube kits are pouched with an ENFit® Single-Port to Y-Port Adapter and Pinch Clamp.

REPLACEMENT Y-PORT ADAPTERS		
DESCRIPTION	ORDER #	BOX/QTY
ENFit® Single-Port to Y-Port Adapter	4-0005-ISOSAF	10


Traditional Non-Balloon G-Tube, Non-Balloon

Adjustable Length Gastrostomy Tubes

Capsule Monarch® G – Tube



Stages of Capsule Bolster Deployment



- 1: Encapsulated Bolster Inserted into Established Stoma
- 2: Tether Pulled to Split Capsule Material, Stiffener Removed

Exclusive A M T Capsule Technology Ease of Placement.
Improve Comfort + Durability. Reduce Leaks.

DESCRIPTION	ENFit® ORDER #
12F	9-C1210-I
14F	9-C1410-I
18F	9-C1810-I
20F	9-C2010-I
12Fw/Guidewire (0.018" x 70 cm)	9-CG1210-I
14Fw/Guidewire (0.018" x 70 cm)	9-CG1410-I
18Fw/Guidewire (0.038" x 90 cm)	9-CG1810-I
20F w/Guidewire (0.038" x 90 cm)	9-CG2010-I

AMT Capsule Monarch® G-Tube is Packaged 1/box.
All ENFit® Capsule Monarch® kits are pouched with an ENFit® Single-Port to Y-Port Adapter and Pinch Clamp. For Guidewire Kits: Guidewire is included in a separate pouch, sterile, length varies by French size

Replacement Y-Port Adapters



ENFit®
Single-Port Adapter

ENFit®
Single-Port to Y-Port
Adapter

Spare adaptors available please see p16

G-JET® Gastric-Jejunal Feeding Tube Family



PAEDIATRIC

micro G-JET®
Low Profile Gastric-Jejunal
Button

ENFit® | 14F Gastric – 8F
Jejunal



LOW PROFILE

G-JET® Button
Low Profile Gastric-Jejunal Button

ENFit® | 14F, 16F, 18F



TRADITIONAL

Traditional G-JET®
Adjustable Length Gastric-Jejunal
Enteral Feeding Tube


ENFit® | 14F, 16F, 18F, 22F

micro G-JET® Button

Paediatric G-J Enteral Feeding Tube

One of a Kind!

Tubing transitions from a 14F Gastric segment to an 8F Jejun segment.




Extension + Feed Sets Included in Kit

‡AMT Jejunal feed sets are compatible with Mic-Key® devices.



Designed for Paediatric Patients⁵
The Only Low Profile Gastric-Jejunal Feeding Device with an 8F Jejunal Segment



GASTRIC 14F	JEJUNAL 8F	ENFIT® ORDER #
0.8 cm	10 cm	MGJ-1408-10-I
	15 cm	MGJ-1408-15-I
	22 cm	MGJ-1408-22-I
1.0 cm	10 cm	MGJ-1410-10-I
	15 cm	MGJ-1410-15-I
	22 cm	MGJ-1410-22-I
1.2 cm	10 cm	MGJ-1412-10-I
	15 cm	MGJ-1412-15-I
	22 cm	MGJ-1412-22-I
1.5 cm	10 cm	MGJ-1415-10-I
	15 cm	MGJ-1415-15-I
	22 cm	MGJ-1415-22-I
1.7 cm	10 cm	MGJ-1417-10-I
	15 cm	MGJ-1417-15-I
	22 cm	MGJ-1417-22-I
2.0 cm	10 cm	MGJ-1420-10-I
	15 cm	MGJ-1420-15-I
	22 cm	MGJ-1420-22-I
2.3 cm	10 cm	MGJ-1423-10-I
	15 cm	MGJ-1423-15-I
	22 cm	MGJ-1423-22-I
2.5 cm	10 cm	MGJ-1425-10-I
	15 cm	MGJ-1425-15-I
	22 cm	MGJ-1425-22-I


micro G-JET® Button comes packaged as 1/box (sterile). | Kit contents include: device (1), gastric extension set (1), jejunal extension set (1), guide-wire adapter (1), 5cc/10cc balloon syringe (1), 10cc contrast syringe (1), 35cc bolus syringe (1), winged luer adapter (1), and gauze.

G-JET® Button


Low Profile Gastric-Jejunal Enteral Feeding Tube

An Advantage in Enteral Nutrition

Anti-Kink. Exclusive AMT Balloon. Glow Green Technology.



Same French Size, Larger Lumens
Designed to Reduce Clogging



GASTRIC Lumen Comparison⁶

AMT 14F is 58% Larger than Competitor 14F

AMT 16F is 41% Larger than Competitor 16F

AMT 18F is 34% Larger than Competitor 18F

JEJUNAL Lumen Comparison⁶


AMT 14F is — 56% Larger than Competitor 14F
56% Larger than Competitor 16F

AMT 16F is — 131% Larger than Competitor 16F
81% Larger than Competitor 18F

AMT 18F is — 172% Larger than Competitor 18F
57% Larger than Competitor 22F

Mutually Exclusive Ports Leading the Way in Patient Safety

Accidentally feeding the stomach instead of the jejunum can be dangerous! Mutually exclusive **GASTRIC** and **JEJUNAL** ports are designed to reduce the risk of misconnections. Different feed/extension sets are required for each port, so it is imperative to maintain adequate stock of both sets. Always remember to **Feed the Green!**



G-JET® Button

Low Profile Gastric-Jejunal Enteral Feeding Tube

G-JET® Button comes packaged as 1/box (sterile). | Kit contents include: device (1), gastric extension set (1), jejunal extension set (1), guidewire adapter (1), 5cc/10cc balloon syringe (1), 10cc contrast syringe (1), 35cc bolus syringe (1), winged luer adapter (1), and gauze.

Extension + Feed Sets Included in Kit

‡AMT Jejunal feed sets are compatible with Mic-Key® devices.




GASTRIC	JEJUNAL	14F	16F	18F
LENGTH (cm)		ENFIT® ORDER #	ENFIT® ORDER #	ENFIT® ORDER #
1	15	GJ-1410-15-I	GJ-1610-15-I	
	22	GJ-1410-22-I	GJ-1610-22-I	
1.2	15	GJ-1412-15-I	GJ-1612-15-I	
	22	GJ-1412-22-I	GJ-1612-22-I	GJ-1812-22-I
	30	GJ-1412-30-I	GJ-1612-30-I	GJ-1812-30-I
1.5	15	GJ-1415-15-I	GJ-1615-15-I	
	22	GJ-1415-22-I	GJ-1615-22-I	GJ-1815-22-I
	30	GJ-1415-30-I	GJ-1615-30-I	GJ-1815-30-I
	45	GJ-1415-45-I	GJ-1615-45-I	GJ-1815-45-I
1.7	15	GJ-1417-15-I	GJ-1617-15-I	
	22	GJ-1417-22-I	GJ-1617-22-I	GJ-1817-22-I
	30	GJ-1417-30-I	GJ-1617-30-I	GJ-1817-30-I
	45	GJ-1417-45-I	GJ-1617-45-I	GJ-1817-45-I
2	22	GJ-1420-22-I	GJ-1620-22-I	GJ-1820-22-I
	30	GJ-1420-30-I	GJ-1620-30-I	GJ-1820-30-I
	45	GJ-1420-45-I	GJ-1620-45-I	GJ-1820-45-I
2.3	22	GJ-1423-22-I		
	30	GJ-1423-30-I	GJ-1623-30-I	GJ-1823-30-I
	45	GJ-1423-45-I	GJ-1623-45-I	GJ-1823-45-I
2.5	22	GJ-1425-22-I		
	30	GJ-1425-30-I	GJ-1625-30-I	GJ-1825-30-I
	45	GJ-1425-45-I	GJ-1625-45-I	GJ-1825-45-I
2.7	30	GJ-1427-30-I	GJ-1627-30-I	GJ-1827-30-I
	45	GJ-1427-45-I	GJ-1627-45-I	GJ-1827-45-I
3	45	GJ-1430-45-I	GJ-1630-45-I	GJ-1830-45-I
3.5	45	GJ-1435-45-I	GJ-1635-45-I	GJ-1835-45-I
4	45	GJ-1440-45-I	GJ-1640-45-I	GJ-1840-45-I
4.4	45	GJ-1444-45-I	GJ-1644-45-I	GJ-1844-45-I
5	45	GJ-1450-45-I	GJ-1650-45-I	GJ-1850-45-I
5.5	45	GJ-1455-45-I	GJ-1655-45-I	GJ-1855-45-I
6	45	GJ-1460-45-I	GJ-1660-45-I	GJ-1860-45-I
6.5	45	GJ-1465-45-I	GJ-1665-45-I	GJ-1865-45-I
7	45	GJ-1470-45-I	GJ-1670-45-I	GJ-1870-45-I

Traditional Length G-JET


Adjustable Length Gastric-Jejunal Enteral Feeding Tube

Comprehensive French Sizes

Designed for both pediatrics and adults ranging from 14F – 22F.



The G – JET® You Trust, Now Offered in Traditional Length Designed to Reduce Clogging



GASTRIC Lumen Comparison⁶

AMT 14F is Exclusive, Not Offered by Competitor

AMT 16F is 41% Larger than Competitor 16F

AMT 18F is 34% Larger than Competitor 18F

JEJUNAL Lumen Comparison⁶

AMT 14F is — 56% Larger than Competitor 16F
— 22% Larger than Competitor 18F

AMT 16F is — 131% Larger than Competitor 16F
— 81% Larger than Competitor 18F

AMT 18F is — 172% Larger than Competitor 18F
— 57% Larger than Competitor 22F


JEJUNAL LENGTH (cm)	GASTRIC 14 F ENFIT® ORDER #	GASTRIC 16F ENFIT® ORDER #	GASTRIC 18F ENFIT® ORDER #	GASTRIC 22F ENFIT® ORDER #
15	TGJ-14-15-I	TGJ-16-15-I		
22	TGJ-14-22-I	TGJ-16-22-I	TGJ-18-22-I	
30	TGJ-14-30-I	TGJ-16-30-I	TGJ-18-30-I	TGJ-22-30-I
45	TGJ-14-45-I	TGJ-16-45-I	TGJ-18-45-I	TGJ-22-45-I

Traditional G-JET® comes packaged as 1/box (sterile). | Kit contents include: device (1), 5cc/10cc balloon syringe (1), 35cc bolus syringe (1), and gauze. NOTE: G-JET® devices should never be cut to custom lengths.

ENFit® Extension Sets


Decompression – Medication

8-XXXX-ISOSAF extension sets are compatible with:
MiniONE® Gastric Port (G-Port) + G-JET® Jejunal Port (J-Port)



	DESCRIPTION	ORDER CODE	BOX
	2" Straight Single ENFit® Medication Set w/ Twist Cap	8-0211-I	10
	12" Straight Single ENFit® Feed Set w/ Twist Cap	8-1211-I	10
	12" Right Angle Single ENFit® Feed Set w/ Twist Cap	8-1221-I	10
	24" Straight Purple Bolus Feed Set	8-2411-I	10
	24" Right Angle Purple Bolus Feed Set	8-2421-I	10
	12" ENFit® Port Right Angle w/ Twist Cap	8-1226-I	10
	12" Right Angle Single ENFit® Feed Set	8-1223-I	10
	24" Right Angle Single ENFit® Feed Set	8-2423-I	10
	2" Straight Dual ENFit® Y-Port Medication Set	8-0212-I	10
	12" Straight Dual ENFit® Y-Port Feed Set	8-1212-I	10
	2" Right Angle Dual ENFit® Y-Port Medication Set	8-0222-I	10
	12" Right Angle Dual ENFit® Y-Port Feed Set	8-1222-I	10
	24" Right Angle Dual ENFit® Y-Port Feed Set	8-2422-I	10
	12" Right Angle Dual ENFit® Y-Port Feed Set	8-1255-I	10
	24" Right Angle Dual ENFit® Y-Port Feed Set	8-2455-I	10
	12" Right Angle With Y-Port	8-1224-I	10
	24" Right Angle With Y-Port	8-2424-I	10
	12" ENFit® Port Right Angle Double Lock Single	8-1223D-I	10
	24" ENFit® Port Right Angle Double Lock Single	8-2423D-I	10

6-XXXX-I extension sets are compatible with:
Mini Classic Buttons & G-JET® Gastric Port (G-Port)



	DESCRIPTION	ORDER CODE	BOX
	2" Straight Single ENFit® Medication Set w/ Twist Cap	6-0211-I	10
	12" Straight Single ENFit® Extension Set w/ Twist Cap	6-1211-I	10
	24" Straight Single ENFit® Extension Set w/ Twist Cap	6-2411-I	10
	12" Right Angle Single ENFit® Extension Set w/ Twist Cap	6-1221-I	10
	24" Right Angle Purple Bolus Extension Set w/ Twist Cap	6-2421-I	10
	12" Right Angle Single ENFit® Extension Set w/ Twist Cap	6-1226-I	10
	12" Right Angle ENFit® Extension Set	6-1223-I	10
	24" Right Angle Single ENFit® Extension Set	6-2423-I	10
	2" Straight Clear Dual ENFit® Y-Port Medication Set	6-0212-I	10
	12" Straight Clear Dual ENFit® Y-Port Medication Set	6-1212-I	10
	12" Right Angle Clear Dual ENFit® Y-Port Extension Set	6-1222-I	10
	24" Right Angle Clear Dual ENFit® Y-Port Extension Set	6-2422-I	10
	12" Right Angle Medication Set	6-1224-I	10
	24" Right Angle Medication Set	6-2424-I	10
	12" ENFit® Port Right Angle Double Lock Single	6-1223D-I	10
	24" ENFit® Port Right Angle Double Lock Single	6-2423D-I	10


CINCH®

Tube Securement Device


From 4F to 30F, The CINCH® Has You Covered

Hypoallergenic


Gentle adhesive for sensitive skin.



MEDIUM CINCH®



LARGE CINCH®




AMT CINCH®




MEDIUM CINCH® TUBE SECUREMENT DEVICE		
DESCRIPTION	ORDER #	BOX/QTY
Medium CINCH® 4F – 17F	CINCH417M-10	10
Medium CINCH® 4F – 17F	CINCH417M-50	50
LARGE CINCH® TUBE SECUREMENT DEVICE		
DESCRIPTION	ORDER #	BOX/QTY
Large CINCH® 9F – 30F	CINCH930L-10	10
Large CINCH® 9F – 30F	CINCH930L-50	50

The AMT CINCH® was designed with the most delicate skin in mind. Our hypoallergenic silicone adhesive provides a gentle yet secure attachment. Secures 4F to 30F medical tubes.

Front



Back



The CINCH® comes fully assembled with an easy-peel split backing.

AMT Accessories



The AMT Clamp™ device can reduce disconnects between a “Christmas Tree”/ Stepped adapter and a feed set.

AMT CLAMP SECUREMENT DEVICE		
DESCRIPTION	ORDER #	BOX/QTY
Tube Securement Device	4-3000	10

I.P. KIT ACCESSORIES - AVAILABLE SEPARATELY	
ORDER #	DESCRIPTION
1804	Balloon Stoma Measuring Device – Guidewire Compatible (up to 0.038”) 
1704	Balloon Stoma Measuring Device – Non-Guidewire Compatible 
4-6000	Introducer/Stiffener – 4.6” Length Guidewire Compatible (up to 0.038”). Device Compatibility: 12F and larger Low Profile Devices 
4-5000	Introducer/Stiffener – Non-Guidewire Compatible. Device Compatibility: 12F and larger MINI Button Devices 

IP Kit accessories are individually pouched and sold as 10 per box (sterile)



ENFit® TRANSITION ADAPTERS		
DESCRIPTION	ORDER #	BOX/QTY
Male ENFit® Transition Adapter	TRN101	10
ENFit® Transition Stepped Adapter	TRN102	10
Female ENFit® Transition Adapter	TRN201	10
ENFit® Bolus Transition Adapter	TRN202	10
ENFit® Y-Port Transition Adapter	TRN203	10
ENFit® Y-Port Adapter w/ Enhanced Med-Port	TRN204	10



ENFit®
Single-Port
Adapter




ENFit®
Single-Port to
Y-Port Adapter

TRADITIONAL LENGTH G-TUBE REPLACEMENT ADAPTERS		
DESCRIPTION	ORDER #	BOX/QTY
Capsule Dome, Capsule Monarch®		
12F ENFit® Single-Port Adapter	4-2012-ISOSAF	10
14F ENFit® Single-Port Adapter	4-2014-ISOSAF	10
20F ENFit® Single-Port Adapter	4-2020-ISOSAF	10
24F ENFit® Single-Port Adapter	4-2024-ISOSAF	10
Balloon, Capsule Dome, Capsule Monarch®		
ENFit® Single-Port to Y-Port Adapter	4-0005-ISOSAF	10

Initial Placement Kits

Initial Placement Gastrostomy Sets

AMT Quality, Right from the Start.



1. Multiple Size Dilators
For use with 10F – 18F devices.



2. Ease of Placement
An AMT exclusive, size-specific, guidewire compatible Introducer/Stiffener.



3. Ensure a Proper Fit
Use AMT's Balloon Stoma Measuring Device with each placement.



Designed to Minimize Trauma Associated with Stoma Tract Development

- Tapered Dilators Help Ensure Proper Dilation
- Requires Only One Endoscopic Pass
- Compatible with Laparoscopic, Endoscopic, and Fluoroscopic Placement
- For Use with 10F – 18F Devices

BASIC I.P. KITS†	
ORDER CODE	KIT CONTENTS
IP-DIL	(4) Dilators (7F, 12F, 16F, 20F)
	(1) Guidewire, 45 cm (0.038")
	(2) Introducer Needles, 18 Gauge (short, long)
	(1) Guidewire Compatible Introducer/Stiffener

I.P. Kits are packaged 1/box (sterile).

COMPREHENSIVE I.P. KITS*	
ORDER CODE	KIT CONTENTS
IP-1218	(5) Dilators (7F, 12F, 16F, 18F, 22F)
	(1) Guidewire, 90 cm (0.038")
	(2) Introducer Needles, 18 Gauge (short, long)
	(1) Guidewire Compatible BSMD
	(2) Guidewire Compatible Introducer/Stiffener
	(1) Safety Scalpel, #11
	(1) 5 cc Syringe
IP-1218-S	(4) 4" x 4" Gauze
ORDER CODE	KIT CONTENTS
IP-1218-S	(5) Dilators (7F, 12F, 16F, 18F-S, 22F)
	(1) Guidewire, 90 cm (0.038")
	(2) Introducer Needles, 18 Gauge (short, long)
	(1) Guidewire Compatible BSMD
	(2) Guidewire Compatible Introducer/Stiffener
	(1) Safety Scalpel, #11
	(1) 5 cc Syringe
	(4) 4" x 4" Gauze

I.P. Kits are packaged 1/box (sterile).

I.P. Kits include one sheathed dilator, noted in bold above.

The IP-1812-S is a peel away sheath

Bladder and Bowel Management

MiniSPC – The first and only licenced Low Profile Suprapubic Catheter

Length (cm)	12F	14F	16F	18F	20F	24F	12F	14F	16F	18F	20F	24F
0.8	SPC-1208-12	SPC1408-12					SPC-1208-24	SPC1408-24				
1	SPC-1210-12	SPC1410-12	SPC-1610-12	SPC-1810-12	SPC-2010-12		SPC-1210-24	SPC1410-24	SPC-1610-24	SPC-1810-24	SPC-2010-24	
1.2	SPC-1212-12	SPC1412-12	SPC-1612-12	SPC-1812-12	SPC-2010-12		SPC-1212-24	SPC1412-24	SPC-1612-24	SPC-1812-24	SPC-2010-24	
1.5	SPC-1215-12	SPC1415-12	SPC-1615-12	SPC-1815-12	SPC-2015-12	SPC-2415-12	SPC-1215-24	SPC1415-24	SPC-1615-24	SPC-1815-24	SPC-2015-24	SPC-2415-24
1.7	SPC-1217-12	SPC1417-12	SPC-1617-12	SPC-1817-12	SPC-2017-12	SPC-2417-12	SPC-1217-24	SPC1417-24	SPC-1617-24	SPC-1817-24	SPC-2017-24	SPC-2417-24
2	SPC-1220-12	SPC1420-12	SPC-1620-12	SPC-1820-12	SPC-2020-12	SPC-2420-12	SPC-1220-24	SPC1420-24	SPC-1620-24	SPC-1820-24	SPC-2020-24	SPC-2420-24
2.3	SPC-1223-12	SPC1423-12	SPC-1623-12	SPC-1823-12	SPC-2023-12	SPC-2423-12	SPC-1223-24	SPC1423-24	SPC-1623-24	SPC-1823-24	SPC-2023-24	SPC-2423-24
2.5	SPC-1225-12	SPC1425-12	SPC-1625-12	SPC-1825-12	SPC-2025-12	SPC-2425-12	SPC-1225-24	SPC1425-24	SPC-1625-24	SPC-1825-24	SPC-2025-24	SPC-2425-24
2.7	SPC-1227-12	SPC1427-12	SPC-1627-12	SPC-1827-12	SPC-2027-12	SPC-2427-12	SPC-1227-24	SPC1427-24	SPC-1627-24	SPC-1827-24	SPC-2027-24	SPC-2427-24
3	SPC-1230-12	SPC1430-12	SPC-1630-12	SPC-1830-12	SPC-2030-12	SPC-2430-12	SPC-1230-24	SPC1430-24	SPC-1630-24	SPC-1830-24	SPC-2030-24	SPC-2430-24
3.5	SPC-1235-12	SPC1435-12	SPC-1635-12	SPC-1835-12	SPC-2035-12	SPC-2435-12	SPC-1235-24	SPC1435-24	SPC-1635-24	SPC-1835-24	SPC-2035-24	SPC-2435-24
4	SPC-1240-12	SPC1440-12	SPC-1640-12	SPC-1840-12	SPC-2040-12	SPC-2440-12	SPC-1240-24	SPC1440-24	SPC-1640-24	SPC-1840-24	SPC-2040-24	SPC-2440-24
4.4		SPC1444-12	SPC-1644-12	SPC-1844-12	SPC-2044-12	SPC-2444-12		SPC1444-24	SPC-1644-24	SPC-1844-24	SPC-2044-24	SPC-2444-24
5		SPC1450-12	SPC-1650-12	SPC-1850-12	SPC-2050-12	SPC-2450-12		SPC1450-24	SPC-1650-24	SPC-1850-24	SPC-2050-24	SPC-2450-24
5.5		SPC1455-12	SPC-1655-12	SPC-1855-12	SPC-2055-12	SPC-2455-12		SPC1455-24	SPC-1655-24	SPC-1855-24	SPC-2055-24	SPC-2455-24
6		SPC1460-12	SPC-1660-12	SPC-1860-12	SPC-2060-12	SPC-2460-12		SPC1460-24	SPC-1660-24	SPC-1860-24	SPC-2060-24	SPC-2460-24
6.5		SPC1465-12						SPC1465-24				
7		SPC1470-12						SPC1470-24				
8		SPC1480-12						SPC1480-24				
9		SPC1490-12						SPC1490-24				
10		SPC14100-12						SPC14100-24				
12" Drainage Set – Box 10 – S1-1221							24" Drainage Set – Box 10 – S1-2421					




MiniSPC™ Button offers an improved quality of life compared to traditional urinary drainage devices.

The detachable drainage sets maximize discretion, allowing for greater mobility and ability to perform daily tasks.


Bladder and Bowel Management MiniACE® Enema Button

Low Profile Antegrade Continence Enema Button

A New Standard of Care in Bowel Management
The Only FDA Cleared Balloon Button Device for ACE Procedures




Irrigation Set Included in Kit



MiniACE® Button comes packaged as 1/box (sterile). | Kit contents include: device (1), irrigation set (1), 5cc syringe (1), introducer/stiffener (1), gauze, and water soluble lubricant (1).

The Difference is Clear:
MiniACE® Enema Button: Low Profile Antegrade Continence Enema Button



	10F	12F	14F
LENGTH (cm)	ENFIT® ORDER #	ENFIT® ORDER #	ENFIT® ORDER #
0.8	ACE-1008-I	ACE-1208-I	ACE-1408-I
1	ACE-1010-I	ACE-1210-I	ACE-1410-I
1.2	ACE-1012-I	ACE-1212-I	ACE-1412-I
1.5	ACE-1015-I	ACE-1215-I	ACE-1415-I
1.7	ACE-1017-I	ACE-1217-I	ACE-1417-I
2	ACE-1020-I	ACE-1220-I	ACE-1420-I
2.3	ACE-1023-I	ACE-1223-I	ACE-1423-I
2.5	ACE-1025-I	ACE-1225-I	ACE-1425-I
2.7	ACE-1027-I	ACE-1227-I	ACE-1427-I
3	ACE-1030-I	ACE-1230-I	ACE-1430-I
3.5	ACE-1035-I	ACE-1235-I	ACE-1435-I
4	ACE-1040-I	ACE-1240-I	ACE-1440-I
4.4	ACE-1044-I	ACE-1244-I	ACE-1444-I
5	ACE-1050-I	ACE-1250-I	ACE-1450-I
5.5	ACE-1055-I	ACE-1255-I	ACE-1455-I
6	ACE-1060-I	ACE-1260-I	ACE-1460-I
6.5	ACE-1065-I	ACE-1265-I	ACE-1465-I
7	ACE-1070-I	ACE-1270-I	ACE-1470-I
8	ACE-1080-I	ACE-1280-I	ACE-1480-I
9	ACE-1090-I	ACE-1290-I	ACE-1490-I
10	ACE-10100-I	ACE-12100-I	ACE-14100-I

Contact Kebomed

*Need Help with a Product or
Have a Question?*

GET IN TOUCH, WE'RE HERE TO HELP!

At Kebomed our team is committed to empowering patients, carers and Healthcare professionals through educational resources and support.

Our trained Clinical Sales Consultants are available to offer meetings/ appointments in person or virtually. Additionally, our Customer Services team are ready to answer your questions.



Phone:
01884 266377



EMAIL:
Enteral.sales@kebomed.co.uk

Procurement and Healthcare professionals can place an order via email.



SET UP AN ACCOUNT:

New customers will need to submit a New Account Application prior to placing their first order.



PLACE YOUR ORDER:

Items that are in stock will be sent out that day if placed before 12 noon or next day if after 12 Noon

Customer Service team is available from 9:00 am to 5:00 pm

References:

1. Seder, C.W., et al. Crit Care Med 2010;38:797-801
2. Badosa, E.L. Nutr Hosp. 2017;34(4):907-913
3. Claim based on internal testing of 5F – 18F devices.
4. Abdelhadi, Ruba A., Dekonenko, C., Dorman, Robert M., et al. A Prospective, Randomized Cross-over Trial of Two Low-profile Gastrostomy Buttons to Determine Family Preference. JPGN 2020;70:386 – 388.
5. The micro G-JET® is cleared for patients > 6 kg.
6. Claim based on internal testing of 14-22F devices.
7. Claim based on internal testing of 16F devices.
8. Hill D.S., O'Neill, J.K., Powell, R.J., et al. Surgical smoke—a health hazard in the operating theatre: a study to quantify exposure and a survey of the use of smoke extractor systems in UK plastic surgery units. Journal of Plastic Reconstructive & Aesthetic Surgery. 2012; 65(7):911-916.
9. Pierce J.S., Lacey S.E., Lippert, J.F., et al. Laser-Generated Air Contaminants from Medical Laser Applications: A State-of-the-Science Review of Exposure Characterization, Health Effects, and Control. JOEH. 2011; 8:447-466
10. Tomita, Y., Mihashi, S., Nagata, K., et al. Mutagenicity of smoke condensates induced by CO2-laser irradiation and electrocauterization. Mutation research 89 2 (1981):145-149.
11. Hill D.S., O'Neill, J.K., Powell, R.J., et al. Surgical smoke—a health hazard in the operating theatre: a study to quantify exposure and a survey of the use of smoke extractor systems in UK plastic surgery units. Journal of Plastic Reconstructive & Aesthetic Surgery. 2012; 65(7):911-916.
12. Based on internal testing of 14F, 18F, 20F, and 24F devices.
13. Claim based on internal testing of 10F devices.