



LiNA CoolBrain

Surgical headband with cooling effect

LiNA CoolBrain

The LiNA CoolBrain is a pre-moistened surgical headband with an automatic cooling effect to keep the surgeon and other OR staff cool and dry during surgery. Furthermore, the CoolBrain is soft and comfortable to wear and equipped with one-size elastic band to fit all head sizes.

The CoolBrain is made of PVA Foam which has an extremely large surface. The large surface allows the moisture to evaporate, giving the desired cooling effect. During surgery, the foam will absorb any perspiration, to continue the cooling.

The CoolBrain comes sterile and packed in a special aluminum pouch retaining moisture inside the headband for more than 3 years from date of production. The moisture consists of 99.8% water and 0.2% preservative.

- *Specially designed aluminium pouch retains moisture in headband even after 3 years*
- *Pre-moistened with one-size elastic band*
- *Ready for use*
- *Sterile*



Ordering information:

LiNA CoolBrain

REF.	Product	Packing
CB-200	LiNA CoolBrain - Surgical Headband made of PVA foam with pore size 130 micron. Pre-moistened and ready for use. Single use.	Box of 25, sterile