



Date: 2 September 2025

Memo: AMT Statement regarding: Feed Set Compatibility – to ensure optimal safety for patients

Dear Healthcare Professional

To ensure optimal performance and reliability, AMT recommends the exclusive use of AMT extension sets with AMT devices.

Although other brands may appear compatible, AMT cannot guarantee their safety, performance, or durability when used in place of AMT extension sets.

AMT's extension sets are designed for safe use with AMT products, not just to be compatible.

Some **non**-AMT sets, especially those with white tips, negate a built-in safety feature. Our Glow Green™ technology was created to protect users from misconnections and to make nighttime feeding easier for parents and carers.

Any questions relating to the above please contact your local Clinical Sales Consultant or our customer services team on 01884 266377, Option 2