





Innovation in Gynaecology




LiNA McCartney Tube™


*Improving the outcome of
Minimally Invasive Gynaecology*




Identifies where to safely divide the uterine vessels



Provides a clear template for safe and easy colpotomy



Maintains pneumoperitoneum and provides an easy exit for specimens




Delivers a convenient entry and exit for needle and suture and presents the vaginal edges for laparoscopic closure of the vaginal vault



LiNA McCartney Tube™

The LiNA McCartney Tube offers you a unique single use device, sterile and ready to use.



safe

- Maintains pneumoperitoneum after colpotomy
- Protection of the ureters from trauma during surgery
- Identifies where to safely seal and divide the uterine vessels
- Provides a clear template for bladder dissection and colpotomy
- Improved visualisation designed for a safer procedure

functional

- **Lid with 5mm & 10mm valves**
Prevents loss of pneumoperitoneum and enables insertion of instruments.
- **Fast and simple laparoscopic suturing technique of the vaginal vault**
Delivers a convenient entry and exit for needle and suture and presents the vaginal edges for laparoscopic closure.
- **Handle to control tube positioning during procedures**



versatile

- **Easy removal of specimen, e.g. the uterus, lymph nodes and ovarian cysts**
- **Tip withstands both electrosurgical and ultrasonic energy**
Soft curvature designed for atraumatic insertion and anatomic fit.



Ordering Information:

LiNA McCartney Tube™

35-TV	LiNA McCartney Tube - 35mm diameter, sterile, single use	Box of 5
45-TV	LiNA McCartney Tube - 45mm diameter, sterile, single use	Box of 5